

version number: V1.0

#### SHENZHEN WELL.D MEDICAL ELECTRONICS CO., LTD

Add: Qing Lan 3rd Rd, Well.d Industrial Park, Biopharmaceutical Base, Pingshan New Area, Shenzhen,

Tel: +86 755 36900018

Email: export@welld.com.cn

For More Information, please visit us at www.welld.com.cn. Shenzhen Well.d Medical Electronics Co.,ltd. All right Reserved All Precautions can be referenced from User Manual.












Back fat thickness / loin eye muscle area



10.1" medical high-resolution LCD with full touch screen

### **∠** LEADING TECHNOLOGY, EXCELLENT IMAGING



Automatic identification of probe interface



Synthetic aperture beamforming technology (SA)



Multi-beam parallel processing technology



Pulse Inversion Harmonic Imaging Technology (PHI)



Speckle Noise Suppression Technology



**Precise Transmit** Focus Technology

### Optional probes



Convex probes

Linear probes







Microconvex probe Animal rectal probe

Back-fat probe






## Light

1.65KG Only Equipped with Leather bag


# **Easy operation**

Internal lithium battery, working time over 3 hours Support wireless images transmission



### + Clear Diagnostic Images


"Mimage" Latest core technology embeded




Back fat



Beef ribeye



Pig loineye



Early pregnancy









**Follicles**